CERASORB[®] CPC

DE	Gebrauchsanweisung	2
GB/EN	Instructions for use	8
DK	Brugsanvisning	14
ES	Instrucciones de uso	20
GR	Οδηγίες χρήσης	26
FI	Käyttöohje	32
FR	Instructions d'utilisation	38
IT	Istruzioni per l'uso	44
NL	Gebruiksaanwijzing	50
NO	Bruksanvisning	56
PL	Instrukcja stosowania	62
PT	Instruções de utilização	68
SE	Bruksanvisning	74

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Instructions for use

Description

CERASORB® CPC is a synthetic, biocompatible, osteoconductive and bioresorbable bone substitute material for filling non-load-bearing bone defects.

The material is applied directly to the bone defect from the syringe or with a cannula without any further preparation.

Composition

CERASORB® CPC is a mineral bone cement composed of synthetic calcium and phosphate salts finely dispersed in a biocompatible oil phase made from short-chain triglycerides (caprylic / capric triglycerides) as well as two emulsifiers (polyoxyl-35-castor oil and cetyl phosphate). Caprylic / capric triglycerides and polyoxyl-35-castor oil are made from plant-based raw materials.

The setting reaction of CERASORB® CPC starts after application upon contact with an aqueous solution, such as bodily fluid.

CERASORB® CPC sets in situ to a microcrystalline, calcium-deficient hydroxyapatite (CDHA) and alpha-tricalcium phosphate, which are the main phase. The minor phase consists of monetite and calcite.

The chemical composition and crystalline structure of the final reaction product (CDHA) correspond largely to the mineral components of natural bone.

CERASORB® CPC does not contain any substances of animal origin, added preservatives or pharmacologically active agents. A volume of 1 ml CERASORB® CPC contains approx. 24 mg potassium in the form of K_2HPO_4 (see dosage and application).

Components + percentage (%)	
alpha- tricalcium phosphate	48.35 - 49.88
calcium hydrogen phosphate (monetite)	20.95 - 21.61
calcium carbonate (calcite)	8.06 - 8.32
tricalcium ortho-phosphate (CDHA)	3.23 - 3.33
dipotassium hydrogen phosphate	2.41 - 2.49
Miglyol 812 (caprylic / capric triglycerides)	11.56 - 13.68
Kolliphor ELP (poly-oxyl-35-castor oil)	2.11 - 2.50
Amphisol A (cetyl phosphate)	0.70 - 0.82

Intended use

CERASORB® CPC is a self-setting, synthetic bone substitute material for filling non-infected bone defects.

Areas of application

CERASORB® CPC is intended for filling or reconstructing non-load-bearing bone defects or for filling bone defects that have been sufficiently stabilised by appropriate means.

Particular areas of application are:

- metaphyseal defect fractures, e.g. fractures of the tibia, radius and humerus
- bone defects following resection of benign tumours and cysts
- bone defects after removal or replacement of osteosynthetic implants
- support for the fixation of osteosynthetic implants (e.g. of bone screws)

CERASORB® CPC must be applied only to a well vascularised and non-infected bone site.

The correct repositioning and stabilisation of fractures is to be assured by means of appropriate fixation.

Use

CERASORB® CPC is an implantable product and designed for use in open or minimally-invasive applications.

The CERASORB® CPC package includes a cannula. When using the cannula, the force required for application is increased. It must be noted that, depending on the size of the cannula, some of the CERASORB® CPC will remain in the cannula and is thus not available for filling the defect. The cannula must be shortened according to the requirements (e.g. with Cooper scissors).

CERASORB® CPC can also be applied without using a cannula.

The physician is responsible for the patient's treatment plan, including the duration and timing of the clinical and radiological follow-up. The patient must adhere to the physician's treatment plan. During the pre-operation discussion, the patient should be informed about the treatment conditions with CERASORB® CPC according to the instructions for use.

Preparation of the 0.5 ml, 1 ml and 3 ml version

The CERASORB® CPC package contains a pre-filled syringe and a cannula. Open the packages that contain the syringe with CERASORB® CPC and remove the blue cap. CERASORB® CPC can then be applied to the prepared bone defect directly or using the cannula. To do so, apply slowly, even pressure to the syringe plunger.

Preparation of version 6 ml and 12 ml

The CERASORB® CPC package contains a pre-filled syringe, a spindle drive and a cannula. Open the packages that contain the syringe with CERASORB® CPC. Then open the package with the spindle drive (spindle nut and threaded spindle). The spindle nut must be pushed onto the back end of the barrel until a clear click sound is heard. Check that both sides of the spindle nut are locked.

Screw the threaded spindle in the spindle nut until it is in contact with the piston. Remove the blue cap from the syringe. CERASORB® CPC can then be applied to the prepared bone defect directly or using the cannula by slowly rotating the spindle.

Setting behaviour / hardening

CERASORB® CPC is set after combination or contact with an aqueous solution, such as bodily fluids. Prior to contact with an aqueous solution, CERASORB® CPC can be processed indefinitely.

Like all calcium phosphate cements, CERASORB® CPC must be inserted only after final repositioning and stabilising of the bone defect to prevent manipulation during the further hardening phase and impairing the setting behaviour.

Once it has been applied to the bone defect, CERASORB® CPC must not be manipulated (no dabbing, no corrective measures, no testing of the hardening, etc.).

Since CERASORB® CPC is hardened with the surrounding fluid, the strength development depends on the shape and size of the filled defect. The surface dimensional stability of CERASORB® CPC is achieved after approx. 10 to 15 minutes in an aqueous environment, after approx. 4 days, the compressive strength reaches up to 35 MPa.

CERASORB® CPC takes several days to set in the implantation site and reaches a compressive strength of up to 45 MPa in a fully hardened state.

Dosage and application

CERASORB® CPC is available in the following quantities:

Content per syringe		Potassium content per syringe
0.5 ml	1 g	approx. 12 mg
1 ml	2 g	approx. 24 mg
3 ml	6 g	approx. 70 mg
6 ml	12 g	approx. 140 mg
12 ml	24 g	approx. 280 mg

The amount of CERASORB® CPC required to fill a defect depends on the defect size.

In order to avoid delays during the surgical procedure, it should be ensured prior to treatment that a sufficient number of CERASORB® CPC packages are available to fill the bone defect completely.

However, there must not be a gap of more than 5 minutes between each application in order to avoid disintegration in the hardening phase between the individual doses or of the material already inserted.

The implanted amount of 42 g (21 ml) CERASORB® CPC must not be exceeded for an adult during one surgery.

For minimally-invasive applications, suitable imaging methods must be used to monitor the filling of the bone defect.

Contraindications

CERASORB® CPC must **not** be used in case of:

- acute or chronic infections at the implantation site, e.g. osteomyelitis
- bone defects due to malignant tumours
- bone defects in the area of open epiphyseal plates
- augmentation in the area of the spine, since there is no clinical experience to date
- cranioplasty
- pregnant or nursing women
- known intolerance to an ingredient of CERASORB® CPC (see composition)

CERASORB® CPC must be used only after carefully weighing the risks and benefits in the case of:

- bone metabolism disorders
- endocrinopathies
- immunosuppressive therapy
- simultaneous treatment with medication that has an effect on bone metabolism

Use during pregnancy or lactation

No clinical data for use in pregnant or lactating women is available. Therefore, CERASORB® CPC should not be used during pregnancy or lactation.

Use in children or geriatric patients

The implanted amount of 6 g (3 ml) CERASORB® CPC must not be exceeded for a child during one surgery.

There are no known reasons that would restrict use in elderly patients.

Side effects

No side effects attributable to the application of CERASORB® CPC are known.

Interactions

If the patient is being treated simultaneously with resorption-inhibiting substances (particularly bisphosphonates, NSAIDs – non-steroidal anti-inflammatory drugs), slower resorption of the implant material should be assumed.

No additional interactions with other medical devices or medicinal products are known, as long as they do not have a direct effect on bone metabolism.

CERASORB® CPC is a non-metallic, non-conducting and non-magnetic bone substitute material and therefore is designated as MRI safe.

Precautions and warnings

The use of CERASORB® CPC is restricted to specialists who are familiar with handling bone substitute materials, the appropriate surgical techniques, and the treatment of bone defects from their training.

When using CERASORB® CPC, as with all surgical interventions, surgery-related risks may occur, particularly rejection reactions, pseudarthrosis, swelling, seroma and haematoma formation, pain, and wound healing disorders. Subsequent revision operations may become necessary.

Patients with weaker immune systems (e.g. those suffering from rheumatism or diabetes), smokers and alcoholics are at a higher risk of infections and implant failure. Such patients must be informed by medical staff of the possible risks before surgery.

The treatment of post-operative infections may be hampered by the presence of an implanted foreign body and it may prove necessary to remove the implanted material.

Conventional surgical tools that are part of the standard equipment in the operating theatre can be used to remove bone cement.

Any implanted foreign bodies can lead to allergic or inflammation reactions or fever.

Very rare cases of allergic reactions and anaphylactic shock have been reported in literature associated with the component polyoxyl 35 castor oil, which is contained in CERASORB® CPC.

For patients with severely impaired kidney function, adrenal insufficiency or liver cirrhosis, the additional intake of small amounts of potassium can increase the danger of hyperkalaemia or aggravate the existing condition. This is also the case for patients with reduced renal potassium excretion induced by the use of medicines (e.g. heparin, ACE inhibitors, potassium-sparing diuretics, spironolactone, non-steroidal anti-inflammatory drugs, cyclosporin A).

Since potassium is only successively released from CERASORB® CPC and the amount of potassium contained in CERASORB® CPC corresponds to only a fraction of the amount present in the daily nutritional intake, a low risk can be assumed, even with severely impaired kidney function. In particular when large amounts of implanted CERASORB® CPC are used, serum potassium levels must be monitored.

For patients whose kidneys function normally, there is no risk even when large quantities (max. 21 ml) of CERASORB® CPC are implanted.

CERASORB® CPC must only be implanted after sufficient debridement of the bone defect in order to ensure a vital bone site. The defect must be completely filled in order to establish direct osseous contact between CERASORB® CPC and the surrounding bone.

When inserting CERASORB® CPC into a bone defect, it must be ensured that the material does not escape into adjacent soft tissue or blood vessels. To prevent an embolism, it must be ensured that no material is inserted into open veins or arteries, particularly when inserting CERASORB® CPC under pressure into defects that are enclosed on all sides. When CERASORB® CPC is inserted into heavily bleeding bone defects, the material may be pushed out again through bleeding pressure. It is recommended to monitor bleeding before applying CERASORB® CPC and to continue filling the defect only when bleeding is under control.

Due to its mechanical properties, CERASORB® CPC can support the stabilisation of bone defects, but the actual stabilisation must be ensured by other means.

CERASORB® CPC must not be mixed with aqueous solutions, including those of autologous or allogenic origin (e.g. blood) as this results in changes to the material properties of CERASORB® CPC.

CERASORB® CPC is resorbed by biological processes and replaced by the body's own bone. CERASORB® CPC resorbs slowly. Depending on the implantation conditions and the metabolic activity at the implantation

site CERASORB® CPC can also remain permanently in the body as an osseous integrated material.

Any unused content of opened or damaged packages must not be used for further operations and must be discarded (see disposal).

Shelf life

The expiry date is on the labels. The product must not be used after this date.

Storage

CERASORB® CPC should be stored at room temperature between 5°C and 25°C. The required application force and setting speed may change if stored outside of the stated temperature range.

Short-term deviations from the storage conditions during transportation (maximum of 5 days) have been tested for the range from -23°C to +55°C and do not result in an unusable product.

As this is a sterile medical device, the packaging must be protected from damage to avoid contamination of the product itself.

Syringes from damaged packages must not be used.

Sterilisation procedure

CERASORB® CPC is a sterile medical device. It is sterilised using gamma radiation. Due to the risk of infection transmission and/or the potential impairment of product performance, CERASORB® CPC must not be cleaned or resterilised.

CERASORB® CPC is intended for single use only.

Disposal

No special disposal measures are required for unopened products. Expired or contaminated material must be disposed of in accordance with the hospital standard disposal procedure.

Information






For further information, please contact your supplier or the manufacturer.

Users and/or patients should report any serious incident in relation to the product to the manufacturer and the competent authority of the Member State in which the user and/or patient is located. The Summary of Safety and Clinical Performance (SSCP) is published on the INNOTERE website and is available at the following LINK: www.innotere.de/downloads.

Responsible manufacturer

INNOTERE GmbH
Meissner Str. 191
01445 Radebeul
Germany
+49 351 2599 9410
www.innotere.de

Symbols

	(Notified Body: TÜV Rheinland LGA Products GmbH)
	Catalogue number
	Batch code
	Quantity
	Use-by date
	Manufacturer
	Consult instructions for use
	Caution, consult accompanying documents
	Sterilized using gamma irradiation
	Do not re-sterilize
	Do not use if package is damaged
	Do not re-use
	Temperature limitation
	Unique Device Identification
	Distributor

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